


Anti-Ptpn1 Polyclonal Antibody For Wb, Elisa - P18031

Item Number: CM0023718



Introduction

Validated Rabbit Polyclonal Antibody targeting PTPN1 (PTP1B). Suitable for WB and ELISA; cross-reacts with Human and Rat. Detects a 50kDa protein. Ideal for studying tyrosine phosphatase signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PTPN1 Rabbit pAb
Remarks/Alias	PTP1B; N1
Species	Human
GeneID (Human)	5770
GeneID	5770
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 281-435 of human PTPN1 (NP_002818.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P18031
Protein Weight	50kDa
Shipping	Ice bag