


Anti-Ptpn2 Polyclonal Antibody (Ko Validated) For Wb, If/Icc, Elisa - P17706

Item Number: CM0022048



Introduction

KO validated rabbit polyclonal antibody against human PTPN2 (TCPTP), suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Involved in T-cell signaling, cytokine regulation, and glucose homeostasis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PTPN2 Rabbit pAb
Remarks/Alias	PTN2; PTPT; TCPTP; TC-PTP; TCELLPTP; N2
Species	Human
GeneID (Human)	5771
GeneID	5771
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-353 of human PTPN2 (NP_536348.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17706
Protein Weight	48kDa
Shipping	Ice bag