


Anti-Ptpn2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P17706

Item Number: CM0019101



Introduction

PTPN2 Rabbit Polyclonal Antibody (TCPTP) validated for Western Blot, Immunofluorescence/Immunocytochemistry, ELISA in human, mouse, rat. Targets key phosphatase in immune and cancer signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PTPN2 Rabbit pAb
Remarks/Alias	PTN2; PTPT; TCPTP; TC-PTP; TCELLPTP; PTPN2
Species	Human
Gene ID (Human)	5771
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-353 of human PTPN2 (NP_536348.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17706
Protein Weight	48kDa
Shipping	Ice bag