


Anti-Ptpn23 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H3S7

Item Number: CM0030278



Introduction

Rabbit polyclonal antibody targeting human PTPN23 (HD-PTP). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for endosomal sorting and ciliogenesis research.

[Learn More](#)

Parameter	Value
Product Name	PTPN23 Rabbit pAb
Remarks/Alias	HDPTP; HD-PTP; NEDBASS; PTP-TD14; PTPN23
Species	Human
GeneID (Human)	25930
GeneID	25930
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H3S7
Protein Weight	179kDa
Shipping	Ice bag