


Anti-Ptpn3 Polyclonal Antibody For Wb And Elisa - P26045

Item Number: CM0026666



Introduction

High-quality rabbit polyclonal antibody against human PTPN3 (PTPH1), validated for Western blot and ELISA. Detects endogenous PTPN3 in human, mouse, and rat samples. Suitable for studying membrane-cytoskeleton dynamics and phosphatase activity.

[Learn More](#)

Parameter	Value
Product Name	PTPN3 Rabbit pAb
Remarks/Alias	PTPH1; PTP-H1; PTPN3
Species	Human
GeneID (Human)	5774
GeneID	5774
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26045
Protein Weight	104kDa
Shipping	Ice bag