


Anti-Ptprb Polyclonal Antibody For Wb, Elisa - P23467

Item Number: CM0028579



Introduction

Rabbit polyclonal antibody against human PTPRB (VE-PTP), validated for WB and ELISA. Critical regulator of blood vessel remodeling and angiogenesis; ideal for endothelial cell biology research.

[Learn More](#)

Parameter	Value
Product Name	PTPRB Rabbit pAb
Remarks/Alias	PTPB; HPTPB; VEPTP; HPTP-BETA; R-PTP-BETA; PTPRB
Species	Human
GeneID (Human)	5787
GeneID	5787
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1173-1260 of human PTPRB (NP_002828.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P23467
Protein Weight	224kDa
Shipping	Ice bag