


Anti-Ptprd Polyclonal Antibody For Wb, Ihc-P, Elisa - P23468

Item Number: CM0028949



Introduction

PTPRD Rabbit pAb validated for Western blot, IHC-P, and ELISA. Detects human and mouse receptor-type tyrosine-protein phosphatase delta, a key regulator of neuronal synaptic differentiation.

[Learn More](#)

Parameter	Value
Product Name	PTPRD Rabbit pAb
Remarks/Alias	HPTP; PTPD; HPTPD; HPTPDELTA; RPTPDELTA; R-PTP-delta; PTPRD
Species	Human
GeneID (Human)	5789
GeneID	5789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 475-574 of human PTPRD (NP_002830.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P23468
Protein Weight	215kDa
Shipping	Ice bag