


Anti-Ptprf Rabbit Polyclonal Antibody For Wb And Elisa - P10586

Item Number: CM0028452



Introduction

Rabbit polyclonal antibody against human PTPRF (LAR), suitable for WB and ELISA. Cross-reacts with human and mouse. Target is a 213 kDa cell adhesion receptor with protein tyrosine phosphatase activity.

[Learn More](#)

Parameter	Value
Product Name	PTPRF Rabbit pAb
Remarks/Alias	LAR; BNAH2; PTPRF
Species	Human
GeneID (Human)	5792
GeneID	5792
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1016-1250 of human PTPRF (NP_002831.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10586
Protein Weight	213kDa
Shipping	Ice bag