


## Anti-Ptprn2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92932

Item Number: CM0014462



### Introduction

Rabbit polyclonal antibody against human PTPRN2 (IA-2beta, phogrin), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for secretion and diabetes research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | PTPRN2 Rabbit pAb                                                                                                  |
| Remarks/Alias    | IAR; ICAAR; PTPRP; IA-2beta; R-PTP-N2; PTPRN2                                                                      |
| Species          | Human                                                                                                              |
| GeneID           | 5799                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-242 of human PTPRN2 (NP_570858.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q92932                                                                                                             |
| Protein Weight   | 111kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |