


Anti-Ptprr Monoclonal Antibody For Wb, Elisa - Q15256

Item Number: CM0002821



Introduction

Rabbit monoclonal antibody targeting human PTPRR (Receptor-type tyrosine-protein phosphatase R), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for MAPK signaling research.

[Learn More](#)

Parameter	Value
Product Name	PTPRR Rabbit mAb
Remarks/Alias	PTPRQ; EC-PTP; PCPTP1; PTP-SL; PTPBR7; PTPRR
Species	Human
Gene ID (Human)	5801
Gene ID	5801
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PTPRR (Q15256)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q15256
Protein Weight	74kDa
Shipping	Ice bag