


Anti-Ptprr Rabbit Polyclonal Antibody For Wb, Elisa - Q15256

Item Number: CM0015116



Introduction

PTPRR rabbit polyclonal antibody validated for WB and ELISA, with cross-reactivity to human and rat. Targets receptor-type tyrosine-protein phosphatase R, a MAPK regulator. Ideal for signaling and phosphatase research.

[Learn More](#)

Parameter	Value
Product Name	PTPRR Rabbit pAb
Remarks/Alias	PTPRQ; EC-PTP; PCPTP1; PTP-SL; PTPBR7; PTPRR
Species	Human
Gene ID (Human)	5801
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-227 of human PTPRR (NP_002840.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q15256
Protein Weight	74kDa
Shipping	Ice bag