

Anti-Ptpu Rabbit Polyclonal Antibody For Wb, Elisa - Q92729

Item Number: CM0031004



Introduction

Rabbit polyclonal antibody targeting human PTPRU (receptor-type tyrosine-protein phosphatase U). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying cell adhesion, proliferation, and epithelial integrity.

[Learn More](#)

Parameter	Value
Product Name	PTPRU Rabbit pAb
Remarks/Alias	FMI; PTP; PCP-2; PTP-J; PTPRO; PTPU2; PTP-PI; PTP-RO; PTPPSI; hPTP-J; R-PTP-U; R-PTP-PSI; PTPRU
Species	Human
GeneID (Human)	10076
GeneID	10076
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PTPRU (NP_005695.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92729
Protein Weight	162kDa
Shipping	Ice bag