

Anti-Ptrh1 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Y79

Item Number: CM0016199

Introduction

PTRH1 Rabbit pAb is a polyclonal antibody raised against recombinant human PTRH1 protein, suitable for WB and ELISA. Targets peptidyl-tRNA hydrolase involved in ribosome-associated quality control.

[Learn More](#)

Parameter	Value
Product Name	PTRH1 Rabbit pAb
Remarks/Alias	PTH1; C9orf115; PTRH1
Species	Human
GeneID (Human)	138428
GeneID	138428
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-214 of human PTRH1 (NP_001002913.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86Y79
Protein Weight	23kDa
Shipping	Ice bag