

Anti-Pts (6-Pyruvoyl Tetrahydrobiopterin Synthase) Polyclonal Antibody For Wb, Elisa - Q03393

Item Number: CM0028868



Introduction

Rabbit polyclonal antibody against human PTS (6-pyruvoyl tetrahydrobiopterin synthase), suitable for WB and ELISA. Reacts with human, mouse, and rat samples. Useful for BH4 biosynthesis and phenylketonuria research.

[Learn More](#)

Parameter	Value
Product Name	PTS Rabbit pAb
Remarks/Alias	PTPS; PTS
Species	Human
Gene ID (Human)	5805
Gene ID	5805
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03393
Protein Weight	16kDa
Shipping	Ice bag