


Anti-Ptx3 Polyclonal Antibody For Wb And Elisa - P26022

Item Number: CM0015684



Introduction

Rabbit polyclonal antibody recognizing human PTX3 (TSG-14/TNFAIP5), validated for western blot and ELISA. Cross-reacts with mouse and rat. Suitable for studies on innate immunity, inflammation, and pentraxin-related biology.

[Learn More](#)

Parameter	Value
Product Name	PTX3 Rabbit pAb
Remarks/Alias	TSG-14; TNFAIP5; PTX3
Species	Human
Gene ID (Human)	5806
Gene ID	5806
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-381 of human PTX3 (NP_002843.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26022
Protein Weight	42kDa
Shipping	Ice bag