


Anti-Pu.1/Spi1 Monoclonal Antibody For Wb, If, Ihc-P, Fc, Elisa - P17947

Item Number: CM0007005



Introduction

Rabbit monoclonal PU.1/SPI1 antibody validated for Western blot, immunofluorescence, immunohistochemistry, flow cytometry, and ELISA. Targets human PU.1 transcription factor (31 kDa, SWISS: P17947) for hematopoietic research.

[Learn More](#)

Parameter	Specification
Product Name	PU.1/SPI1 Rabbit mAb
Remarks/Alias	OF; PU.1; AGM10; SFPI1; SPI-1; SPI-A; PU.1/SPI1
Species	Human
Gene ID (Human)	6688
Gene ID	6688
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PU.1/SPI1 (NP_003111.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, FC (intra), ELISA
Cross Reactivity	Human
SWISS	P17947
Protein Weight	31kDa
Shipping	Ice bag