


Anti-Puf60 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Uhx1

Item Number: CM0030752



Introduction

Rabbit polyclonal antibody targeting human PUF60, validated for Western blot, immunoprecipitation, and ELISA with cross-reactivity in human, mouse, and rat samples. Ideal for transcription and splicing research.

[Learn More](#)

Parameter	Value
Product Name	PUF60 Rabbit pAb
Remarks/Alias	FIR; VRJS; RoBPI; SIAHBP1; PUF60
Species	Human
GeneID	22827 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 243-542 of human PUF60 (NP_055096.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHX1
Protein Weight	60kDa
Shipping	Ice bag