

Anti-Pumilio 1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q14671

Item Number: CM0003430



Introduction

CamelBio offers the Anti-Pumilio 1 Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human and mouse Pumilio 1 (Q14671). Ideal for RNA-binding protein research. Ships on ice.

[Learn More](#)

Parameter	Value
Product Name	Pumilio 1 Rabbit mAb
Remarks/Alias	PUMH; HSPUM; PUMH1; PUML1; SCA47; Pumilio 1
Species	Human
Gene ID	9698
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 104-203 of human Pumilio 1 (Q14671)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14671
Protein Weight	126kDa
Shipping	Ice bag