

Anti-Pumilio 1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14671

Item Number: CM0025661



Introduction

Rabbit polyclonal antibody against human Pumilio 1 (PUM1), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for post-transcriptional regulation and RNA-binding protein studies.

[Learn More](#)

Parameter	Value
Product Name	Pumilio 1 Rabbit pAb
Remarks/Alias	PUMH; HSPUM; PUMH1; PUML1; SCA47; Pumilio 1
Species	Human
GeneID (Human)	9698
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human Pumilio 1 (NP_055491.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14671
Protein Weight	126kDa
Shipping	Ice bag