


Anti-Pumilio 2 Monoclonal Antibody For Wb, Elisa - Q8Tb72

Item Number: CM0003053



Introduction

Rabbit monoclonal antibody against human Pumilio 2 (PUM2), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying RNA-binding, post-transcriptional regulation, and stem cell biology.

[Learn More](#)

Parameter	Value
Product Name	Pumilio 2 Rabbit mAb
Remarks/Alias	PUMH2; PUML2; Pumilio 2
Species	Human
GeneID (human)	23369
GeneID	23369
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Pumilio 2 (Q8TB72).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TB72
Protein Weight	114kDa
Shipping	Ice bag