

Anti-Purinergic Receptor P2Y6 Monoclonal Antibody For Wb And Elisa - Q15077

Item Number: CM0012249

Introduction

Rabbit monoclonal antibody targeting human P2Y purinoceptor 6 (P2Y6, Q15077). Validated for WB and ELISA applications. Recognizes a 36 kDa GPCR that mediates UDP-induced calcium signaling.

[Learn More](#)

Parameter	Value
Product Name	Purinergic Receptor P2Y6 Rabbit mAb
Remarks/Alias	P2Y6; Purinergic Receptor P2Y6
Species	Human
Gene ID (Human)	5031
Gene ID	5031
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 229-328 of human Purinergic Receptor P2Y6 (Q15077).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15077
Protein Weight	36kDa
Shipping	Ice bag