

Anti-Purinergic Receptor P2Y6 Rabbit Polyclonal Antibody For Wb, Elisa - Q15077

Item Number: CM0018972

Introduction

Rabbit polyclonal antibody against human P2RY6 (P2Y purinoceptor 6). Suitable for WB and ELISA. Cross-reacts with human and mouse. Target is a G protein-coupled receptor for UDP, UTP, and ATP, mediating calcium signaling. Ideal for GPCR and nucleotide signaling research.

[Learn More](#)

Parameter	Value
Product Name	Purinergic Receptor P2Y6 Rabbit pAb
Remarks/Alias	P2Y6; Purinergic Receptor P2Y6
Species	Human
GeneID (Human)	5031
GeneID	5031
Immunogen	Synthetic peptide corresponding to amino acids 200-300 of human P2RY6 (NP_789766.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15077
Protein Weight	36kDa
Shipping	Ice bag