


Anti-Pus1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y606

Item Number: CM0029340



Introduction

High-quality rabbit polyclonal antibody against human PUS1, validated for WB and ELISA. Ideal for studying pseudouridylation, mRNA splicing, and mitochondrial disorders. Cross-reacts with human PUS1 at 47 kDa.

[Learn More](#)

Parameter	Value
Product Name	PUS1 Rabbit pAb
Remarks/Alias	MLASA1; PUS1
Species	Human
Gene ID	80324 (Human)
Immunogen	Recombinant fusion protein corresponding to amino acids 308-427 of human PUS1 (NP_079491.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y606
Protein Weight	47kDa
Shipping	Ice bag