

Anti-Pus10 Polyclonal Antibody For Wb And Elisa - Q3Mit2

Item Number: CM0021548



Introduction

Rabbit polyclonal antibody targeting human PUS10, validated for Western blotting and ELISA. Cross-reacts with mouse and rat. Ideal for investigating tRNA pseudouridylation and miRNA processing pathways.

[Learn More](#)

Parameter	Value
Product Name	PUS10 Rabbit pAb
Remarks/Alias	DOBI; Hup10; CCDC139; PUS10
Species	Human
GeneID (Human)	150962
GeneID	150962
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q3MIT2
Protein Weight	60kDa
Shipping	Ice bag