


Anti-Pus7 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Pz0

Item Number: CM0018937



Introduction

Anti-PUS7 rabbit polyclonal antibody validated for WB and ELISA, targeting human PUS7 with cross-reactivity to mouse and rat. Ideal for studying RNA pseudouridylation and splicing regulation.

[Learn More](#)

Parameter	Value
Product Name	PUS7 Rabbit pAb
Remarks/Alias	IDDABS; PUS7
Species	Human
Gene ID (Human)	54517
Gene ID	54517
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 352-661 of human PUS7 (NP_061915.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PZ0
Protein Weight	75kDa
Shipping	Ice bag