

Anti-Pvrig Rabbit Polyclonal Antibody For Wb, Elisa - Q6Dki7

Item Number: CM0023192



Introduction

High-quality rabbit polyclonal anti-PVRIG (CD112R) antibody, validated for WB and ELISA. Ideal for exploring coinhibitory T-cell signaling and immune checkpoint pathways.

[Learn More](#)

Parameter	Value
Product Name	PVRIG Rabbit pAb
Remarks/Alias	CD112R; C7orf15; PVRIG
Species	Human
GeneID (Human)	79037
GeneID	79037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 193-326 of human PVRIG (NP_076975.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6DKI7
Protein Weight	34kDa
Shipping	Ice bag