

Anti-Pvrl1 Rabbit Polyclonal Antibody For Wb, Elisa - Q15223

Item Number: CM0022530



Introduction

Rabbit polyclonal antibody against human PVRL1 (Nectin-1/HveC/CD111), validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for cell adhesion and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	PVRL1 Rabbit pAb
Remarks/Alias	ED4; PRR; HlgR; HV1S; HVEC; OFC7; PRR1; PVRR; CD111; PVRL1; PVRR1; SK-12; CLPED1; nectin-1
Species	Human
GeneID (Human)	5818
GeneID	5818
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-517 of human PVRL1 (NP_002846.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15223
Protein Weight	57kDa
Shipping	Ice bag