

Anti-Pvrl1/Nectin1 Monoclonal Antibody For Fc - Q15223

Item Number: CM0007721



Introduction

Rabbit monoclonal antibody targeting human Nectin-1 (PVRL1/CD111), validated for flow cytometry. Suitable for investigating cell adhesion, herpesvirus entry, and neural development. Reacts with human.

[Learn More](#)

Parameter	Specification
Product Name	PVRL1/NECTIN1 Rabbit mAb
Remarks/Alias	ED4; PRR; HlgR; HV1S; HVEC; OFC7; PRR1; PVRR; CD111; PVRL1; PVRR1; SK-12; CLPED1; nectin-1; PVRL1/NECTIN1
Species	Human
GeneID (human)	5818
GeneID	5818
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human PVRL1/NECTIN1 (NP_002846.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q15223
Protein Weight	57kDa/51kDa/39kDa
Shipping	Ice bag