


Anti-Pvrl3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nqs3

Item Number: CM0029931



Introduction

Polyclonal antibody raised against human PVRL3 (Nectin-3, CD113). Suitable for Western blot and ELISA, cross-reacts with mouse. Useful for studies of cell adhesion, axon guidance, and cell junctions. - Q9NQS3

[Learn More](#)

Parameter	Value
Product Name	PVRL3 Rabbit pAb
Remarks/Alias	PPR3; PRR3; CD113; PVRL3; PVRR3; CDW113; NECTIN-3
Species	Human
GeneID (Human)	25945
GeneID	25945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 127-366 of human PVRL3 (NP_001230215.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NQS3
Protein Weight	61kDa
Shipping	Ice bag