

Anti-Pvrl4 (Nectin-4) Rabbit Polyclonal Antibody For Wb And Elisa - Q96Ny8

Item Number: CM0019303



Introduction

Rabbit polyclonal antibody against human PVRL4 (Nectin-4, LNIR, PRR4) validated for WB and ELISA. Suitable for cell adhesion and measles virus research. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	PVRL4 Rabbit pAb
Remarks/Alias	LNIR; PRR4; EDSS1; PVRL4; nectin-4
Species	Human
Gene ID (Human)	81607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-146 of human PVRL4 (NP_112178.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96NY8
Protein Weight	55kDa
Shipping	Ice bag