

Anti-Pwp1 Monoclonal Antibody For Wb, Elisa - Q13610

Item Number: CM0003358

Introduction

This rabbit monoclonal antibody specifically recognizes human PWP1 (Periodic tryptophan protein 1 homolog), a nucleolar chromatin-associated factor. Validated for WB and ELISA; recombinant immunogen corresponds to aa 460-501. Predicted MW 56 kDa. UniProt Q13610.

[Learn More](#)

Parameter	Value
Product Name	PWP1 Rabbit mAb
Remarks/Alias	NCLB; IEF-SSP-9502; PWP1
Species	Human
Gene ID (Human)	11137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-501 of human PWP1 (NP_008993.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13610
Protein Weight	56kDa
Shipping	Ice bag