

Anti-Pxdn Rabbit Monoclonal Antibody For Wb, Elisa - Q92626

Item Number: CM0012187

Introduction

Rabbit monoclonal antibody against human PXDN for WB and ELISA applications. Cross-reacts with mouse and rat. Ideal for extracellular matrix, basement membrane, and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	PXDN Rabbit mAb
Remarks/Alias	PXN; VPO; MG50; PRG2; ASGD7; COPOA; D2S448; D2S448E; hsPxd01
Species	Human
GeneID (Human)	7837
GeneID	7837
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 431-520 of human PXDN (NP_036425.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92626
Protein Weight	165kDa
Shipping	Ice bag