

Anti-Pxdn Rabbit Polyclonal Antibody For Wb And Elisa - Q92626

Item Number: CM0022472

Introduction

Rabbit polyclonal antibody targeting human PXDN (Peroxidasin homolog). Validated for WB and ELISA, cross-reacts with mouse. Suitable for extracellular matrix and basement membrane research.

[Learn More](#)

Parameter	Value
Product Name	PXDN Rabbit pAb
Aliases	PXN, VPO, MG50, PRG2, ASGD7, COPOA, D2S448, D2S448E, hsPxd01, PXDN
Species	Human
Gene ID	7837
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1320-1479 of human PXDN (NP_036425.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92626
Protein Weight	165 kDa
Shipping	Ice bag