

Anti-Pxk Rabbit Polyclonal Antibody For Wb And Elisa - Q7Z7A4

Item Number: CM0014064

Introduction

PXK Rabbit Polyclonal Antibody targets the PX domain-containing protein kinase-like protein (MONAKA), a regulator of Na,K-ATPase and synaptic transmission. Suitable for WB and ELISA in human, mouse, and rat samples. Detects a 65kDa band.

[Learn More](#)

Parameter	Value
Product Name	PXK Rabbit pAb
Remarks/Alias	SLOB; MONAKA; PXK
Species	Human
Gene ID (Human)	54899
Gene ID	54899
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PXK (NP_060241.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z7A4
Protein Weight	65kDa
Shipping	Ice bag