


Anti-Pxr Polyclonal Antibody For Wb, Elisa - 075469

Item Number: CM0014755



Introduction

Rabbit polyclonal antibody targeting human PXR (NR1I2), validated for WB and ELISA. Detects a 50 kDa protein; suitable for xenobiotic metabolism and drug transport research. Cross-reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	PXR Rabbit pAb
Remarks/Alias	BXR; PAR; PRR; PXR; SAR; SXR; ONR1; PAR1; PAR2; PARq
Species	Human
GeneID (human)	8856
GeneID	8856
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human PXR (NP_071285.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75469
Protein Weight	50 kDa
Shipping	Ice bag