

Anti-Pycr1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P32322

Item Number: CM0008024



Introduction

High-quality rabbit monoclonal antibody against human PYCR1, validated for WB, IF, IHC, and ELISA. Ideal for studying proline metabolism, mitochondrial function, and oxidative stress response. Backed by CamelBio's one-stop IVD raw material supply chain.

[Learn More](#)

Parameter	Value
Product Name	PYCR1 Rabbit mAb
Remarks/Alias	P5C; P5CR; PRO3; PYCR; PIG45; PP222; ARCL2B; ARCL3B
Species	Human
GeneID (Human)	5831
GeneID	5831
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 211-310 of human PYCR1 (NP_008838.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P32322
Protein Weight	33kDa
Shipping	Ice bag