

Anti-Pycr1 Rabbit Polyclonal Antibody For Wb, Elisa - P32322

Item Number: CM0017363



Introduction

High-quality anti-PYCR1 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat PYCR1. Ideal for proline biosynthesis and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	PYCR1 Rabbit pAb
Remarks/Alias	P5C; P5CR; PRO3; PYCR; PIG45; PP222; ARCL2B; ARCL3B; PYCR1
Species	Human
Gene ID (Human)	5831
Gene ID	5831
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-319 of human PYCR1 (NP_008838.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32322
Protein Weight	33kDa
Shipping	Ice bag