


Anti-Pycr2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96C36

Item Number: CM0018190



Introduction

Rabbit polyclonal anti-PYCR2 antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects a 34 kDa protein and is suitable for studying proline biosynthesis and oxidative stress in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	PYCR2 Rabbit pAb
Remarks/Alias	HLD10; P5CR2; PYCR2
Species	Human
GeneID (Human)	29920
GeneID	29920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 271-320 of human PYCR2 (NP_037460.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96C36
Protein Weight	34kDa
Shipping	Ice bag