

## Anti-Pycrl Polyclonal Antibody For Wb And Elisa - Q53H96

Item Number: CM0022882

### Introduction

PYCRL rabbit polyclonal antibody validated for WB and ELISA, with cross-reactivity to mouse. Ideal for investigating proline biosynthesis and metabolic pathways. Optimized for human samples.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | PYCRL Rabbit pAb                                                                                                   |
| Remarks/Alias    | PYCRL                                                                                                              |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 65263                                                                                                              |
| GeneID           | 65263                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 13-286 of human PYCRL (NP_075566.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Mouse                                                                                                              |
| SWISS            | Q53H96                                                                                                             |
| Protein Weight   | 29kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |