

Anti-Pygl Rabbit Polyclonal Antibody - P06737

Item Number: CM0024931



Introduction

Rabbit polyclonal antibody targeting human PYGL (liver glycogen phosphorylase) for WB, IHC-P, IF/ICC, IP, and ELISA applications. Detects human, mouse, and rat PYGL. Ideal for glucose homeostasis and glycogen metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	PYGL Rabbit pAb
Remarks/Alias	GSD6; PYGL
Species	Human
Gene ID (Human)	5836
Gene ID	5836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 690-847 of human PYGL (NP_002854.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06737
Protein Weight	97kDa
Shipping	Ice bag