

Anti-Pygm Rabbit Polyclonal Antibody For Wb, Elisa - P11217

Item Number: CM0028537



Introduction

Rabbit polyclonal antibody targeting human PYGM (myophosphorylase), validated for WB and ELISA. Reacts with human, mouse, and rat samples. For glycogen metabolism and McArdle disease research.

[Learn More](#)

Parameter	Value
Product Name	PYGM Rabbit pAb
Remarks/Alias	GSD5; PYGM
Species	Human
Gene ID	5837
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 743-842 of human PYGM (NP_005600.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11217
Protein Weight	97kDa
Shipping	Ice bag