

Anti-Pym1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Brp8

Item Number: CM0022726

Introduction

High-quality anti-PYM1 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat PYM1 (Q9BRP8), a key exon junction complex regulator involved in NMD and translation.

[Learn More](#)

Parameter	Value
Product Name	PYM1 Rabbit pAb
Remarks/Alias	PYM; WIBG; PYM1
Species	Human
GeneID (human)	84305
GeneID	84305
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-204 of human PYM1 (NP_115721.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRP8
Protein Weight	23kDa
Shipping	Ice bag