


Anti-Pyrophosphatase 1 Monoclonal Antibody For Wb And Elisa - Q15181

Item Number: CM0003512



Introduction

High-specificity rabbit monoclonal antibody targeting human Pyrophosphatase 1 (PPA1, Q15181), suitable for Western blotting and ELISA. Cross-reacts with mouse and rat, enabling multi-species metabolic and signaling research.

[Learn More](#)

Parameter	Value
Product Name	Pyrophosphatase 1 Rabbit Monoclonal Antibody
Remarks/Alias	PP, PP1, IOPPP, HEL-S-66p, SID6-8061, Pyrophosphatase 1
Species	Human
GeneID (Human)	5464
GeneID	5464
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Pyrophosphatase 1 (Q15181).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15181
Protein Weight	33 kDa
Shipping	Ice bag