


Anti-Pyruvate Carboxylase (Pc) Rabbit Mab For Wb And Elisa - P11498

Item Number: CM0012782



Introduction

Rabbit monoclonal antibody against human Pyruvate carboxylase (PC), validated for WB and ELISA. Detects a 130 kDa band in human, mouse, and rat samples. Suitable for metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Pyruvate carboxylase (PC) Rabbit mAb
Remarks/Alias	PCB; PC; Pyruvate carboxylase (PC)
Species	Human
GeneID (Human)	5091
GeneID	5091
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Pyruvate carboxylase (PC) (P11498)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P11498
Protein Weight	130kDa
Shipping	Ice bag