


Anti-Pyruvate Carboxylase (Pc) Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P11498

Item Number: CM0026265



Introduction

Rabbit polyclonal antibody against human Pyruvate carboxylase (PC). Validated for WB, IF/ICC, IP, and ELISA applications. Cross-reacts with mouse and rat PC. Suitable for gluconeogenesis and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Pyruvate carboxylase (PC) Rabbit pAb
Remarks/Alias	PCB; Pyruvate carboxylase (PC)
Species	Human
Gene ID (Human)	5091
Gene ID	5091
Immunogen	Recombinant fusion protein containing amino acids 879-1178 of human Pyruvate carboxylase (NP_001035806.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11498
Protein Weight	130kDa
Shipping	Ice bag