


Anti-Qars Rabbit Polyclonal Antibody For Wb, Ihc, If, Elisa - P47897

Item Number: CM0025347



Introduction

Rabbit polyclonal antibody targeting human QARS (GlutaminyI-tRNA synthetase). Validated for WB, IHC, IF, ELISA. Reacts with human and mouse. Suitable for brain development and protein synthesis research.

[Learn More](#)

Parameter	Value
Product Name	QARS Rabbit pAb
Remarks/Alias	QARS; GLNRS; MSCCA; PRO2195
Species	Human
GeneID (Human)	5859
GeneID	5859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human QARS (NP_005042.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P47897
Protein Weight	88kDa
Shipping	Ice bag