

Anti-Qdpr Polyclonal Antibody For Wb And Elisa - P09417

Item Number: CM0028189



Introduction

Rabbit polyclonal antibody against human QDPR (DHPR, SDR33C1), validated for WB and ELISA. Reacts with human, mouse, rat. Suitable for tetrahydrobiopterin biosynthesis and phenylketonuria research.

[Learn More](#)

Parameter	Value
Product Name	QDPR Rabbit pAb
Remarks/Alias	DHPR; PKU2; HDHPR; SDR33C1; QDPR
Species	Human
GeneID (human)	5860
GeneID	5860
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-244 of human QDPR (NP_000311.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09417
Protein Weight	26kDa
Shipping	Ice bag