


Anti-Qki Monoclonal Antibody For Wb - Q96Pu8

Item Number: CM0006292



Introduction

QKI Rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Targets the KH domain-containing RNA-binding protein QKI, a key regulator of RNA splicing, stability, and myelination.

[Learn More](#)

Parameter	Value
Product Name	QKI Rabbit mAb
Remarks/Alias	QK; Hqk; QK1; QK3; hqkl; QKI
Species	Human
Gene ID (human)	9444
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PU8
Protein Weight	38kDa
Shipping	Ice bag