


## Anti-Qki Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q96Pu8

Item Number: CM0002472



### Introduction

QKI Rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Recognizes human, mouse, rat QKI protein (~38kDa). Ideal for RNA metabolism and myelination research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | QKI Rabbit mAb                                                                                                      |
| Remarks/Alias    | QK; Hqk; QK1; QK3; hqkl; QKI                                                                                        |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 9444                                                                                                                |
| GeneID           | 9444                                                                                                                |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 242-341 of human QKI (Q96PU8). |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, IF/ICC, IP, ELISA                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q96PU8                                                                                                              |
| Protein Weight   | 38kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |