


Anti-Qpct Polyclonal Antibody For Wb, If/Icc, Elisa - Q16769

Item Number: CM0025233



Introduction

Rabbit polyclonal antibody against human QPCT for Western blot, IF/ICC, and ELISA.

Recombinant immunogen; cross-reacts with mouse, rat. Ideal for studying pyroglutamyl peptides and amyloid-beta.

[Learn More](#)

Parameter	Value
Product Name	QPCT Rabbit pAb
Remarks/Alias	QC; GCT; sQC; QPCT
Species	Human
GeneID (human)	25797
GeneID	25797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-361 of human QPCT (NP_036545.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16769
Protein Weight	41kDa
Shipping	Ice bag